

SCH4U Lab: pH of Salts and K_a

Aaron Huang

April 28, 2026

Purpose

1. To measure the pH of salt solutions
2. To determine the K_a of weak acids using pH measurements
3. To prepare and test a buffer solution

Introduction

This lab explores the pH of salt solutions, the acid dissociation of weak acids, and the properties of buffer solutions. When a neutralization occurs, a salt is formed. Neutral salts form when a strong acid and a strong base react. A non-neutral salt forms when a weak acid or a weak base is involved, it exhibits either acidic or basic properties due to hydrolysis reactions in water.

Weak acids/bases only partially dissociate in aqueous solution. The strength of a weak acid or bases is quantified by its acid dissociation constant, K_a or K_b respectively. By measuring the pH of buffer solutions containing a weak acid/weak base and its conjugate base/acid, we can determine the K_a/K_b value. A buffer solution resists changes in pH when small amounts of acid or base are added, making it useful in many chemical and biological applications.

This lab consists of three parts. First, the pH of different salt solutions will be measured using a pH probe. Second, buffer solutions will be prepared by combining a weak acid with its conjugate base, and their pH will be recorded to calculate K_a values. Finally, the buffer's ability to resist pH changes will be tested by adding a strong base and comparing the pH change to that of pure water treated the same way.

Safety considerations: Wear gloves at all times during this experiment, as some solutions may cause skin irritation upon direct contact.

Procedure

Part 1: pH of Salt Solutions

<u>Materials/Equipment</u>	<u>Actions Taken</u>
Beaker (100 mL)	Clean 100 mL beaker was prepared
Graduated Cylinder	20.0 mL of 0.100 M NaCl solution was measured and added to beaker
pH Probe	pH probe was calibrated according to manufacturer instructions pH of NaCl solution was measured and recorded
Deionized Water	pH probe was rinsed thoroughly with deionized water Procedure was repeated for 0.100 M Na ₂ CO ₃ solution Procedure was repeated for 0.100 M NaHSO ₃ solution pH probe was rinsed after each measurement

Part 2: K_a of a Weak Acid

<u>Materials/Equipment</u>	<u>Actions Taken</u>
Beaker (250 mL)	Clean 250 mL beaker was prepared
Graduated Cylinder	20.0 mL of 0.1 M NaHSO_3 solution was measured and added to beaker
	20.0 mL of 0.100 M Na_2SO_3 solution was measured and added to beaker
Glass Rod	Solution was stirred thoroughly to mix
pH Probe	pH of buffer solution was measured and recorded
Deionized Water	pH probe was rinsed thoroughly
	Procedure was repeated using CH_3COOH and NaCH_3COO solutions
	pH of second buffer solution was measured and recorded

Part 3: Buffer Solution Testing

<u>Materials/Equipment</u>	<u>Actions Taken</u>
Buffer Solution (from Part 2)	40.0 mL of $\text{NaHSO}_3/\text{Na}_2\text{SO}_3$ buffer solution from Part 2 was prepared
Graduated Cylinder	10.0 mL of NaOH solution was measured and added to buffer
pH Probe	Resulting pH was measured and recorded
Beaker	40 mL of deionized water was measured and added to clean beaker
Graduated Cylinder	10.0 mL of NaOH solution was measured and added to water
pH Probe	Resulting pH was measured and recorded

Observations

Part 1: pH of Salt Solutions

Table 1: pH Measurements of Salt Solutions

Solution	pH
0.100 M NaCl	
0.100 M Na_2CO_3	
0.100 M NaHSO_3	

Part 2: K_a of Weak Acids

Part 3: Buffer Capacity

Table 2: pH of Buffer Solutions

Buffer	Acid	pH
$\text{NaHSO}_3/\text{Na}_2\text{SO}_3$	HSO_3^-	
$\text{NaCH}_3\text{COO}/\text{CH}_3\text{COOH}$	CH_3COOH	

Table 3: pH Changes Upon Addition of NaOH

Solution	Initial pH	After 10 mL NaOH	pH Change
Buffer ($\text{HSO}_3^-/\text{SO}_3^{2-}$)			
Deionized Water			